

FLEX-300 Series Nitrile Exam Gloves

Features and Benefits

• Nitrile

• Exam Grade: FDA 510K & CE Rated

• Powder Free

• Non-Sterile

• AQL 1.5

• Fentanyl Tested

Chemotherapy Tested

• Low Dermatitis Potential

• ASTM D6319, D6978

• EN 455 Parts 1, 2, 3

Color: Dark BlueTexture: FingertipWeight: 3.0 gramLength: 240mm



Usage









Laboratory





Responders

Specifications







Industrial







Ordering Information ITEM# SIZE **GLOVES PER BOX BOXES PER CASE GLOVES PER CASE** N30LDS Small 100 10 1000 N30LDM Medium 100 10 1000 N30LDL Large 100 10 1000 N30LDXL X-Large 100 10 1000

Note: 200 & 300 count boxes available on request. MOQ applies.

Email: contact@valiantshield.info



FLEX-300 Series Nitrile Exam Gloves

Specifications

OVERVIEW			
Material	Nitrile		
Powder Content	Powder-Free		
Freedom From Holes	1.5 AQL		
Surface	Fingertip Textured		
Weight (Based on Size M)	3.0 gram		
Color	Dark Blue		
Cuff Length	240mm, Standard		
Not Made With Natural Rubber Latex	Yes		
CE registration	Yes		
Chemo Rated (Yes/No)	Yes, in accordance with FDA ASTM D6978		
Fentanyl Rated (Yes/No)	Yes, in accordance with FDA ASTM D6978		
Food Grade Rated (Yes/No)	Yes, per CFR 177.2600		

THICKNESS MEASUREMENTS		
Finger	Spec ≥ 0.08mm +/-2 Typical Value 0.08mm	
Palm	Spec ≥ 0.07mm +/-2 Typical Value 0.06mm	

PHYSICAL PROPERTIES		
Criteria	Before Aging	After Aging
Force at Break	Spec ≥6N Typical Value 7.68N	Spec <u>></u> 6N Typical Value 8.86N

DIMENSIONS				
Item#	Size	Palm Width (mm)	Length	
N30LDS	Small	85 (+/- 5)	240mm / 9.5 inches	
N30LDM	Medium	95 (+/- 5)	240mm / 9.5 inches	
N30LDL	Large	106 (+/- 5)	240mm / 9.5 inches	
N30LDXL	Extra-Large	115 (+/- 5)	240mm / 9.5 inches	



FLEX-300 Series Nitrile Exam Gloves

Chemotherapy Drugs	Concentration	Breakthrough Detection Time in Minutes
Cammstine (BCNU)	3.3mg/ml	11.1
Cisplatin	1.0mg/ml	> 240
Cyclophosphamide (Cytoxan)	20.0mg/ml	> 240
Dacarbazine (DTJC)	10.0mg/ml	> 240
Doxorubicin Hydrochloride	2.0rng/ml	> 240
Etoposide (Toposar)	20.0mg/ml	> 240
Fluorouracil	50.0mg/rnl	> 240
Paclitaxel (Taxol)	6.0mg/ml	> 240
Thiotepa	10.0mg/ml	40.1
Amethopterin Hydrate	25.0mg/ml	> 240
Methotrexate	25.0mg/ml	> 240
Mitomycin C	0.5mg/ml	> 240
Vincristine Sulfate	1.0mg/ml	> 240

^{*}Please Note that the following drugs have low permeation times

Cammstine 11.1 minutes, Thiotepa 40.1 minutes

Fentanyl Permeation as per ASTM D6978				
TEST DRUG AND CONCENTRATION	MINIMUM BREAKTHROUGH DETECTION TIME (Specimen 1/2/3) (minutes)			
Fentanyl Citrate Injection 100mcg/2mL	No Breakthrough up to 240 minutes			